



## BM-02D Sensor pad – Spare part for BM-02

If you want to use the Nanny in multiple beds, use extra sensor pads (item number BM-02D). Install them in each bed and only change the location of the control unit.

### Description

Up to two sensor pads can be used in each bed with a single control unit if the child has grown bigger.

**BM-02D is a medical device.**

### Technical specifications

**dimensions**

305 x 500 x 15 mm

**weight**

1000 g

**material**

PVC-P

**operating conditions**

+5 °C to +35 °C at 30% to 75% rel. humidity

**transport and storage**

0 °C to +40 °C, rel. humidity 10% to 85%